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B|BRAUN
SHARING EXPERTISE

Aesculap Sterile Technology

Check all components of the sterile container visually for damage and correct function before every use:

1. Metal parts not deformed

- Check container base for dents, especially rim of base
- Aluminium lid not warped

2. Check filter retainers

- Filter retainers must be prestressed and have full surface contact at the edge.

- Locating pin for filter retention plate in lid must not be loose

- Filter retainer snaps audibly into position

3. Intactness of seals

- Seals are present and intact
- No cracks, fractures etc.

4. Flawless plastic parts

- No cracks or visible damage to plastic parts, no loose parts

5. Filter inspection

- Filters undamaged (no kinks, holes, racks or crevices)

- Remove cover for inspection (remove anticlockwise, replace clockwise)
- Primeline filters undamaged (no kinks, holes, cracks or crevices)
- Leave filter in cartridge during inspection


- Ensure that the triangles are aligned together (locked position)

6. Locking mechanism fully functional

- The locking mechanism must be fully functional and show no damage
- Lubricate the locking mechanism hinges with Steriliit® from time to time

7. Correct closure of lid

- Locking mechanism of lid snaps audibly into place on the base counterpart

8. Carrying handles undamaged

- Carrying handles are intact and show no visible damage

9. Flawless sterile containers

- Only use sterile containers in flawless condition
- If components are damaged, replace with original spare parts or repair immediately

10. Labelling and sealing

- Sterile container closed
- Either sealed with production label, or production label and plastic seal attached

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

Major article

Sterility maintenance study: Dynamic evaluation of sterilized rigid containers and wrapped instrument trays to prevent bacterial ingress

Harry L. Shaffer MS^{a,*}, Delbert A. Harnish MS^b, Michael McDonald MS^b,
Reid A. Vernon BS^c, Brian K. Heimbuch MS^b

^a Sterilization Consulting Services, LLC, Highlands Ranch, CO

^b Engineering Science Division, Applied Research Associates, Panama City, FL

^c United States Air Force Academy, Colorado Springs, CO

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